


Aesculap Surgical Technologies

Carbotrac



Carbo *trac*

B | BRAUN
SHARING EXPERTISE




Retractors
made from a carbon fibre reinforced polymer

Wundhaken
aus kohlefaserverstärktem Polymer



Carbotrac
Retractor
Wundhaken


BT545P



100 x 40 mm

Carbotrac
Retractor
Wundhaken

BT542P



60 x 40 mm

Carbotrac
Retractor
Wundhaken

BT540P

Retractors

- made from a carbon fibre reinforced polymer
- X-ray translucent
- Non-glare finish
- Autoclaveable 134 °C/2 bar
- High stability
- Feather light – up to 65% weight reduction

Wundhaken

- aus kohlefaserverstärktem Polymer
- Röntgendurchlässig
- Blendfrei
- Autoklavierbar 134 °C/2 bar
- Hohe Stabilität
- Federleicht – bis zu 65% Gewichtsersparnis



AESCULAP®

All it takes to operate.

B | BRAUN
SHARING EXPERTISE

Aesculap AG & Co. KG

Am Aesculap-Platz
78532 Tuttlingen

Phone +49 7461 95-0

Fax +49 7461 95-2600

www.aesculap.de

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogues and pricelists and to take legal actions.